

Kendall 7325

Kendall 7325

Product Specifications

Dimensions

6"H x 10.75"W x 4.5"D

Weight

9.3 pounds

Power Requirements

115 VAC, 50/60 Hz

Compression Type

Sequential, Gradient Pressure

Compression Cycle

11 seconds compression, decompression time based upon Vascular Refill Detection Measurement

Set Pressure

45 mmHg

Sleeve Cooling

Yes

Cycle Monitor

Yes

Bed Hook

Yes

Alarms

Low Pressure, High Pressure, Internal Electronics Malfunction

Kendall Express 9525

Kendall Express 9525

Product Specifications

Dimensions

6.2" x 7"W x 4.5"D

Weight

3.5 pounds

Power Requirements

115 VAC, 50/60 Hz

Compression Type

Sequential, Gradient Pressure

Set Pressure

45 mmHg ankle, 40mmHg calf, 30 mmHg thigh

Cycle Duration

71 variable

Bed Hook

Yes

Alarms

Low Pressure, High Pressure, Internal Electronics Malfunction