

# Devilbiss: Vacu-Aide Portable Suction Pump

## Devilbiss: Vacu-Aide Portable Suction Pump

### Product Specifications

- 9" x 7" x 8"
  - Weighs 6.3 lbs
  - Up to 550 mm Hg vacuum and 27 lpm flow
  - Built-in regulator and vacuum gauge
  - Universal AC power
- 

# Gomco/Allied Medical: Constant/Intermittent Pump Model 6036

## Gomco/Allied Medical: Constant/Intermittent Pump Model 6036

### Product Specifications

#### Vacuum Range

250 mm Hg (constant) / 200 mm Hg (intermittent)

#### Flow Rate

2.5 lpm (constant) / 3.4 (intermittent)

## **Collection Bottle**

Provided with two glass collection bottle systems, one for each mode selection (one 2800 ml and 600 ml)

## **Shipping Weight**

78 lb

## **Dimensions**

**33"H x 23"W x 17"D**

## **Compressor**

### **Diaphragm**

### **Motor**

**1/8 hp Continuous Duty**

## **Features**

- Mode Vacuum Range Adjustable On/Off Cycle Time
- Constant 0-250 mm HG On — 5-90 seconds
- **Intermittent 0-200 mm Hg Off- 23 seconds**

---

# **Gomco/Allied Medical: Portable Aspirator Model 300**

**Gomco/Allied Medical: Portable Aspirator Model 300**

## **Product Specifications**

### **Weight**

23 lbs

**Dimensions**

12"H x 13"W x 9"D

**Compressor**

Diaphragm pump that is designed for general suction use

**Motor**

1/8 Horsepower; creates up to 22" Hg vacuum

---

# **Gomco/Allied Medical: Thoracic Drainage Pump 6050 110/50**

**Gomco/Allied Medical: Thoracic Drainage Pump 6050 110/50**

**Product Specifications**

**Vacuum Range**

0-50 cm H2O

**Flow Rate**

110 LPM minimum

**Dimensions**

33"H x 18"W x 17"D

**Compressor**

## Motor

110/50 Hz

---

# Medela Vario 18 Suction Pump

## Medela Vario 18 Suction Pump

- Vario 18 is both portable and silent and can therefore be used at the patient`s bedside.
- Vario 18 can be used with either a Disposable or a Reusable Collection System.
- Vario 18 offers a range of features that facilitate handling.
- QuatroFlex technology with silent function: can be used in environments in which quiet operation is necessary.
- Precise vacuum settings and fast vacuum build-up: ready when needed.
- One device for multiple needs.
- New users need little instruction to operate.
- Includes limited 2 year warranty

Does not include canister or tubing

*Get the user manual here*

# **Precision Medical: EasyVac Aspirator**

## **Precision Medical: EasyVac Aspirator**

### **Product Specifications**

#### **Dimensions**

Length 13 1/4" (33.7 cm), Width 8" (20.3 cm), Height 9 3/4" (24.8 cm)

#### **Weight**

12.6 lbs. (5.7 kg)

#### **Shipping Weight**

14.5 lbs. (6.6 kg)

#### **Motor**

120 VAC, 60Hz, 5 Amps

#### **Vacuum Range**

0-26 inHg

#### **Horsepower**

1/5

## **Operating Temperature Range**

50°F to 104°F (10°C to 40°C)

## **Storage Conditions**

Temperature: -4°F to 104°F (-20°C to 40°C) Humidity: Max 95% non-condensing

---

# **Schuco: Vac 130 Aspirator**

## **Schuco: Vac 130 Aspirator**

## **Product Specifications**

### **Shipping Weight**

18 pounds

### **Motor**

- 1/10HP, 115V 60Hz

- **Corrosion resistant, thermal protected**

---

# Schuco: Vac 330 Aspirator Suction Unit

## **Schuco: Vac 330 Aspirator Suction Unit**

The Schuco-Vac 330 is a durable medical suction unit which provides variable controlled vacuum regulation between 0-22 Hg and a suction flow rate of 28 lpm. Features vibration-free, easy to read regulator gauge and a hydrophobic/bacteria filter to prevent pump flooding and contamination. Accepts 800cc & 1200cc canisters. 1/10HP motor, 115V 60Hz.

## **Product Specifications**

### **Vacuum Regulation**

0-22 Hg

### **Suction Flow Rate**

28 lpm

### **Collection Canister**

Accepts 800cc and 1200cc

### **Motor**

1/10HP, 115V 60Hz

---

# Schuco: Vac 430 Compact Suction Pump

**Schuco: Vac 430 Compact Suction Pump**

## **Product Specifications**

### **Vacuum Regulation**

0-22 Hg

### **Suction Flow Rate**

28 lpm

### **Collection Canister**

800cc

### **Motor**

1/10HP, 115V 60Hz