

Devilbiss: Vacu-Aide Portable Suction Pump

Devilbiss: Vacu-Aide Portable Suction Pump

Product Specifications

- 9" x 7" x 8"
 - Weighs 6.3 lbs
 - Up to 550 mm Hg vacuum and 27 lpm flow
 - Built-in regulator and vacuum gauge
 - Universal AC power
-

Gomco/Allied Medical: Constant/Intermittent Pump Model 6036

Gomco/Allied Medical: Constant/Intermittent Pump Model 6036

Product Specifications

Vacuum Range

250 mm Hg (constant) / 200 mm Hg (intermittent)

Flow Rate

2.5 lpm (constant) / 3.4 (intermittent)

Collection Bottle

Provided with two glass collection bottle systems, one for each mode selection (one 2800 ml and 600 ml)

Shipping Weight

78 lb

Dimensions

33"H x 23"W x 17"D

Compressor

Diaphragm

Motor

1/8 hp Continuous Duty

Features

- Mode Vacuum Range Adjustable On/Off Cycle Time
 - Constant 0-250 mm HG On — 5-90 seconds
 - **Intermittent 0-200 mm Hg Off- 23 seconds**
-

Gomco/Allied Medical: Portable Aspirator Model 300

Gomco/Allied Medical: Portable Aspirator Model 300

Product Specifications

Weight

23 lbs

Dimensions

12"H x 13"W x 9"D

Compressor

Diaphragm pump that is designed for general suction use

Motor

1/8 Horsepower; creates up to 22" Hg vacuum

Gomco/Allied Medical: Thoracic Drainage Pump 6050 110/50

Gomco/Allied Medical: Thoracic Drainage Pump 6050 110/50

Product Specifications

Vacuum Range

0-50 cm H2O

Flow Rate

110 LPM minimum

Dimensions

33"H x 18"W x 17"D

Compressor

Motor

110/50 Hz

Medela Vario 18 Suction Pump

Medela Vario 18 Suction Pump

- Vario 18 is both portable and silent and can therefore be used at the patient`s bedside.
- Vario 18 can be used with either a Disposable or a Reusable Collection System.
- Vario 18 offers a range of features that facilitate handling.
- QuatroFlex technology with silent function: can be used in environments in which quiet operation is necessary.
- Precise vacuum settings and fast vacuum build-up: ready when needed.
- One device for multiple needs.
- New users need little instruction to operate.
- Includes limited 2 year warranty

Does not include canister or tubing

Get the user manual here

Precision Medical: EasyVac Aspirator

Precision Medical: EasyVac Aspirator

Product Specifications

Dimensions

Length 13 1/4" (33.7 cm), Width 8" (20.3 cm), Height 9 3/4" (24.8 cm)

Weight

12.6 lbs. (5.7 kg)

Shipping Weight

14.5 lbs. (6.6 kg)

Motor

120 VAC, 60Hz, 5 Amps

Vacuum Range

0-26 inHg

Horsepower

1/5

Operating Temperature Range

50°F to 104°F (10°C to 40°C)

Storage Conditions

Temperature: -4°F to 104°F (-20°C to 40°C) Humidity: Max 95% non-condensing

Schuco: Vac 130 Aspirator

Schuco: Vac 130 Aspirator

Product Specifications

Shipping Weight

18 pounds

Motor

- 1/10HP, 115V 60Hz
 - **Corrosion resistant, thermal protected**
-

Schuco: Vac 330 Aspirator Suction Unit

Schuco: Vac 330 Aspirator Suction Unit

The Schuco-Vac 330 is a durable medical suction unit which provides variable controlled vacuum regulation between 0-22 Hg and a suction flow rate of 28 lpm. Features vibration-free, easy to read regulator gauge and a hydrophobic/bacteria filter to prevent pump flooding and contamination. Accepts 800cc & 1200cc canisters. 1/10HP motor, 115V 60Hz.

Product Specifications

Vacuum Regulation

0-22 Hg

Suction Flow Rate

28 lpm

Collection Canister

Accepts 800cc and 1200cc

Motor

1/10HP, 115V 60Hz

Schuco: Vac 430 Compact Suction Pump

Schuco: Vac 430 Compact Suction Pump

Product Specifications

Vacuum Regulation

0-22 Hg

Suction Flow Rate

28 lpm

Collection Canister

800cc

Motor

1/10HP, 115V 60Hz