

Oxygen Flowmeters (FM197A-15L-00)

Please contact us for additional adapter options not shown here.

Pipe Plugs PG-4 or PG-2

1/8 NPT Male Hex Plug

1/4 NPT Male Hex Plug

Puritan Y-Blocks

Air Puritan Y-Block: Male Quick Connect Inlet X Female Coupler Outlets (YO-PQ16)

O2 Puritan Y-Block: Male Quick Connect Inlet X Female Coupler Outlets (YO-PQ24)

VAC/SUC Puritan Y-Block: Male Quick Connect Inlet X Female Coupler Outlets (YO-PQ22)

WAGD/AGSS Puritan Y-Block: Male Quick Connect Inlet X Female Coupler Outlets (YO-PQEV)

Schrader Y-Blocks

Air Schrader Y-Block: Non-swivel Male Quick Connect Inlet X Female Coupler Outlets (YO-SC16)



N2 Schrader Y-Block: Non-swivel Male Quick Connect Inlet X Female Coupler Outlets (YO-SC12)

N20 Schrader Y-Block: Non-swivel Male Quick Connect Inlet X Female Coupler Outlets (YO-SC04)

VAC/SUC Schrader Y-Block: Non-swivel Male Quick Connect Inlet X Female Coupler Outlets (YO-SC22)

O2 Schrader Y-Block: Non-swivel Male Quick Connect Inlet X Female Coupler Outlets (YO-SC24)

Street Tees

1/8 NPT Street Tee All Ports

1/4 NPT Street Tee All Ports

Carbon Dioxide

All hoses are sold by the foot.

The Inner core is made of medical grade non-conductive PVC and is DEHP-Free.

Braid reinforced polyester is also in the structure of the hose.

This hose is only offered in 1/4" Conductive.

Outer cover: Grey Colored PVC

The working pressure: 200 psi max. @ 70° F

Medical Air

All hoses are sold by the foot

The Inner core is made of medical grade non-conductive PVC and is DEHP-Free.

Braid reinforced polyester is also in the structure of the hose.

Outer cover: Yellow Colored PVC

The working pressure: 200 psi max. @ 70° F

The working temperature: 14° F to 150° F

Nitrogen

All hoses are sold by the foot

The Inner core is made of medical grade non-conductive PVC and is DEHP-Free.

Braid reinforced polyester is also in the structure of the hose.

Outer cover: Black Colored PVC

The working pressure: 200 psi max. @ 70° F

The working temperature: 14° F to 150° F

Nitrous Oxide

All hoses are sold by the foot

The Inner core is made of medical grade non-conductive PVC and is DEHP-Free.

Braid reinforced polyester is also in the structure of the hose.

Outer cover: Blue Colored PVC

The working pressure: 200 psi max. @ 70° F

The working temperature: 14° F to 150° F

Oxygen

All hoses are sold by the foot

The Inner core is made of medical grade non-conductive PVC and is DEHP-Free.

Braid reinforced polyester is also in the structure of the hose.

Outer cover: Green Colored PVC

The working pressure: 200 psi max. @ 70° F

The working temperature: 14° F to 150° F

Vacuum(Suction)

All hoses are sold by the foot

The Inner core is made of medical grade non-conductive PVC and is DEHP-Free.

Braid reinforced polyester is also in the structure of the hose.

Outer cover: White Colored PVC

The working pressure: 200 psi max. @ 70° F

The working temperature: 14° F to 150° F

WAGD EVAC

All hoses are sold by the foot

The Inner core is made of medical grade non-conductive PVC and is DEHP-Free.

Braid reinforced polyester is also in the structure of the hose.

Outer cover: Purple Colored PVC

The working pressure: 200 psi max. @ 70° F

The working temperature: 14° F to 150° F