

1043570 Detachable Battery Pack, Philips Trilogy

1043570 Detachable Battery Pack, Philips Trilogy

A high-powered, detachable external battery pack for the Respironics Trilogy ventilator. These batteries are easy to detach and swap. This battery does not require an external charging station. The part number is 1043570.

1055806 Internal Battery Pack, Philips Trilogy

1055806 Internal Battery Pack, Philips Trilogy

The internal battery pack for the Respironics Trilogy ventilator does not require an external charging station. This part should be replaced by a Trained professional. Specifications: Li-Ion 14.8V 4440mAh. The Part Number is 1055806 (1043572).

Abbott/Hospira Plum A+ Infusion Pump Battery

Abbott/Hosira Plum A Infusion Pump Battery

Replacement battery for the Plum A

6 Volt

4.5 Ah

Abbott/Hospira Plum XL Infusion Pump Battery

Abbott/Hospira Plum XL Infusion Pump Battery

Replacement battery for the Plum XL. Also fits the Plum 5000, Plum XLM, & PCA III Infusion Pumps

8 Volt

3.2 Ah

Baxter 6201 Infusion Pump Battery

Baxter 6201 Infusion Pump Battery

Replacement battery for the Baxter 6201

12 Volt

2.0 Ah

Baxter 6301 Infusion Pump Battery

Baxter 6301 Infusion Pump Battery

Replacement battery for the Baxter 6301

12 Volt

3.4 Ah

Baxter AS40/AS50 Syringe Pump Battery

Baxter AS40/AS50 Syringe Pump Battery

Replacement battery for the Baxter AS40/AS50

6 Volt

1.2 Ah

Kangaroo ePump Enteral Feeding Pump Battery

Kangaroo ePump Enteral Feeding Pump Battery

Replacement battery for the Kangaroo ePump

4.8 Volt

3.8 Ah

Masimo Rad-8 Pulse Oximeter Battery

Masimo Rad-8 Pulse Oximeter Battery

Replacement battery for the Rad-8

6 Volt

1.3 Ah

Medfusion 2010 Syringe Pump Battery

Medfusion 2010 Syringe Pump Battery

Replacement battery for the Medfusion 2010. Also fits the 2001 and 2010i

7.2 Volt

1.8 Ah

Medfusion 3500 Syringe Pump Battery

Medfusion 3500 Syringe Pump Battery

Replacement battery for the Medfusion 3500. Also fits the 2500, 3010, 3010a, and Protege

7.2 Volt

2.1 Ah

Nellcor N-595 Pulse Oximeter

Battery

Nellcor N-595 Pulse Oximeter Battery

Replacement battery for the Nellcor N-595. Also fits the NPB Series 190, 195, 290, 295, 390, 395, and 600x

6 Volt

4.5 Ah