

Triumph IP12 Infusion Pump

- Large 4.3 inch water resistant color touch screen
- Dual Battery for providing over 9 hours of battery life
- Enhanced safety features including full screen color coded alarms
- Rate Mode, Time Mode, Body-Weight Mode, Drug Library Mode, and Drip Mode
- Fully customizable drug library with password protection
- Folding carry handle providing easier transport and space saving stackable design
- Universal pole clamp enabling cage mounting
- 10 preinstalled IV lines to include Baxter, Covetrus A Series, Covetrus B Series, ICU Medical, Pivotal A Series, Pivotal B Series, VetOne A Series, VetOne B Series, Zoetis A Series, and Zoetis B Series.

Product Downloads:

- IP12 Infusion Pump Spec Sheet
- IP12 Infusion Pump Quick Start Guide
- IP12 Infusion Pump User Manual

Basic Infusion Startup Process - Triumph IP12 Infusion Pump

Infusion Modes Training Video for Triumph IP12 Infusion Pump and Triumph SP12 Syringe Pump

How To Customize Settings for the Triumph IP12 Infusion Pump and Triumph SP12 Syringe Pump

Sequence Mode Programming for the Triumph IP12 Infusion Pump and Triumph SP12 Syringe Pump

Common User Errors for the Triumph IP12 Infusion Pump

Triumph SP12 Syringe Pump

- Large 4.3 inch water resistant color touch screen
- Dual Battery for providing over 9 hours of battery life
- Enhanced safety features including full screen color coded alarms
- Rate Mode, Time Mode, Body-Weight Mode, Drug Library Mode, and Drip Mode
- Fully customizable drug library with password protection
- Folding carry handle providing easier transport and space saving stackable design
- Universal pole clamp enabling cage mounting
- 6 preinstalled Syringes to include Monoject, BD, Terumo, Nipro, Reli, and Exel

Product Downloads:

- SP12 Syringe Pump Spec Sheet
- SP12 Syringe Pump Quick Start Guide
- SP12 Syringe Pump User Manual

- Acceptable and Current Syringe Brands and Sizes For The Triumph SP12*

*Using a syringe brand and/or size not listed or currently calibrated in the Triumph SP12 can void your warranty.

Basic Infusion Startup Process - Triumph SP12 Syringe Pump

Infusion Modes Training Video for Triumph SP12 Syringe Pump and Triumph IP12 Infusion Pump

How To Customize Settings for the Triumph SP12 Syringe Pump and Triumph IP12 Infusion Pump

How to Correctly Load a Syringe in the Triumph SP12 Syringe Pump

Sequence Mode Programming for the Triumph IP12 Infusion Pump and Triumph SP12 Syringe Pump

Triumph VIP Infusion Pump

- Large 4.3 inch water resistant color touch screen

- **Dual Battery for providing over 9 hours of battery life**
- **Enhanced safety features including full screen color coded alarms**
- **Rate Mode, Time Mode, Body-Weight Mode, Drug Library Mode, and Drip Mode**
- **Fully customizable drug library with password protection**
- **Folding carry handle providing easier transport and space saving stackable design**
- **Universal pole clamp enabling cage mounting**
- **10 preinstalled IV lines to include Baxter, Covetrus A Series, Covetrus B Series, ICU Medical, Pivetal A Series, Pivetal B Series, VetOne A Series, VetOne B Series, Zoetis A Series, and Zoetis B Series.**

Product Downloads:

- [VIP Infusion Pump Spec Sheet](#)
- [VIP Infusion Pump Quick Start Guide](#)
- [VIP Infusion Pump User Manual](#)

Triumph VIP Infusion Pump Basic Setup

Triumph VIP Infusion Pump Features & Settings

Triumph VIP Infusion Pump Drug Library

Triumph VIP Infusion Pump Features During Infusion



Triumph VIP Infusion Pump Extended Warranty Offer

Triumph VIP Pump Common User Error with Power Cord