

Smiths Medical CADD 100 mL Medication Cassette Reservoir 21-7302-24

Smiths Medical CADD 100 mL Medication Cassette Reservoir

CADD Medication Cassette Reservoir, 100 mL (Clear), with Flow Stop, Clamp, Female Leur, Latex Free

All medication cassette reservoirs with flow stop free-flow protection meet the JCAHO National Patient Safety Goals for prevention of free-flow of medication. Flow stop free-flow protection is an integral set-based free-flow protection device that has been added to medication cassette reservoirs.

Order quantity: 12/box

Shelf life: 2 years

Priming volume: 0.2 mL

Rates up to: 125 mL/hr

Tubing length: 8"

Plasticizer: DEHP

Please see the *Instructions for Use* for a complete listing of the indications, contraindications, warnings and precautions.

PRODUCT BENEFITS

- Small
- Lightweight
- Portable

Smiths Medical CADD 250 mL Medication Reservoir

Smiths Medical CADD 250 mL Medication Reservoir

CADD Medication Cassette Reservoir, 250 mL (Clear), with Flow Stop, Clamp, Female Leur, 8" tubing

The CADD™ 250 mL medication cassette reservoirs are designed to provide safer, more secure, ambulatory medication delivery for infusion therapies requiring volumes up to 250 mL. CADD™ medication cassette reservoirs promote patient ambulation and minimize the impact of infusion therapy on a patient's normal lifestyle by decreasing the length of post-operative hospitalizations and time away from home.

Order quantity: 12/box

Shelf life: 2 years

Priming volume: 0.2 mL

Rates up to: 125 mL/hr

Tubing length: 8"

Plasticizer: TOTM

Please see the *Instructions for Use* for a complete listing of the indications, contraindications, warnings and precautions.

PRODUCT BENEFITS

- Small
- Lightweight
- Portable

- Lockbox alternative

Smiths Medical CADD 50 mL Medication Cassette Reservoir 21-7301-24

Smiths Medical CADD 50 mL Medication Cassette Reservoir

CADD Medication Cassette Reservoir, 50 mL (Clear), with Flow Stop, Clamp, Female Leur, Latex Free

All medication cassette reservoirs with flow stop free-flow protection meet the JCAHO National Patient Safety Goals for prevention of free-flow of medication. Flow stop free-flow protection is an integral set-based free-flow protection device that has been added to medication cassette reservoirs.

Order quantity: 12/box

Shelf life: 2 years

Priming volume: 0.2 mL

Rates up to: 125 mL/hr

Tubing length: 8"

Plasticizer: DEHP

Please see the *Instructions for Use* for a complete listing of the indications, contraindications, warnings and precautions.

PRODUCT BENEFITS

- Small
 - Lightweight
 - Portable
-

Smiths Medical CADD Administration Set, 60", Spike

Smiths Medical CADD Administration Set, 60", Spike

Smiths Medical CADD Administration Set, 60" Tubing with Attached Bag Spike, Male Luer, Add-on Anti-siphon Valve, and Clamp

Order quantity: 12/box

Smiths Medical CADD Extension Sets

Smiths Medical CADD Extension Sets

Smiths Medical CADD® Extension Set with male luer, clamp, and integral anti-siphon valve with male luer

Order quantity: 50/box

Priming volume: 1.7 mL

Rates up to: 125 mL/hr

Output tubing length: 60"

Plasticizer: TOTM (non-DEHP)