

Baxter AS50 Pole Clamp

**** Currently Out Of Stock ****

Baxter AS50 Pole Clamp

Pole clamp for the Baxter AS50 Syringe Pump

BiliTx Clear Rubber Sleeve

BiliTx Clear Rubber Sleeve

Clear rubber sleeve for the Respironics BiliTx Phototherapy System.

This piece helps reinforce the connection between the fiber optic panel and the unit.

Philips Part #: 989805622471

(Alternate Part #:1056626)

Disposable Covers for the BiliTx Neonatal Panel

Disposable Covers for the BiliTx Neonatal panel

Disposable covers for the BiliTx Neonatal (short) panel. (Fits Respironics BiliTx Phototherapy Light)

- Provides added comfort during phototherapy treatment and helps prevent cross-contamination between patients
 - 100/box
 - Part Number 121492
 - Phototherapy light not included
-

Disposable Covers for the BiliTx Standard Panel

Disposable Covers for the BiliTx Standard Panel

Disposable covers for the BiliTx Standard (long) panel. (Fits Respironics BiliTx Phototherapy Light)

- Provides added comfort during phototherapy treatment and helps prevent cross-contamination between patients
- 100/box
- Part Number 121493

Maxtec SmartStack IV Pole

Maxtec SmartStack IV Pole

Superior support and durability along with its unique design will allow the SmartStack to stack neatly with other bases of its kind. It is equipped with large, smooth rolling casters offering easy transport in multiple environments.

Product Specifications

- Unique stacking design saves up to 2/3 more space
- Large, smooth rolling caster with 2 wheel locks
- Includes 4 fixed Ram-horn hooks
- Easy, one-handed height adjustment
- Collapsed height of 65"/Extended height of 88"
- 100 lbs max load capacity
- 48 month warranty on squeeze and release mechanism
- Lifetime warranty on pole base

Drop ships directly from manufacturer within 1-7 business days. This product ships at a low flat rate per unit — please purchase it individually, without adding other products to your cart, to get the discount drop ship rate.

Medfusion 3500 Fixed Pole Clamp

Medfusion 3500 Fixed Pole Clamp

Pole clamp for the Medfusion 3500 Syringe Pump

Out Of Stock

Nellcor OEM DOC-10 Patient Cable

Nellcor OEM DOC-10 Patient Cable

OEM Original Nellcor Patient Cable

Features:

- 10 ft. in length
 - Compatible with Nellcor pulse oximeters with OxiMax technology
-

Patient Cable for Smart Monitor 2 Apnea Monitor

Patient Cable for Smart Monitor 2 Apnea

Monitor

- Metal Rydel Connector by Lemo
 - For use with the Smart Monitor 2
-

Respironics Universal Rolling Stand

Respironics Universal Rolling Stand for the V60 Ventilator or BiPAP Vision (Pictured with BiPAP Vision)

This is a Respironics Universal Rolling Stand that can be used for the *V60 Ventilator or BiPAP Vision. There is a different mounting plate for each. Please specify which one you will need when ordering.

**Compatible with V60 ventilators with serial numbers beginning with "1": all serial numbers below 100111553; V60 ventilators with serial numbers beginning with "2": all serial numbers below 201007805.*

Product Specifications

- Built-in humidifier bracket
- O2 hose mount
- Power cord mount
- Accessory basket
- Large castors with brakes

Drop ships directly from manufacturer within 1-7 business days.