

Alaris/Carefusion/Cardinal Medley PCA 8120

Alaris/Carefusion/Cardinal Medley PCA 8120

Product Specifications

Dimensions

4.5"W x 15.0"H x 7.5"D exclusive of security door

Weight

5.5 lbs.

Operating Principle

Positive Displacement

Delivery Units

mg, mcg, mL, mg/h, mcg/h, mL/h

Flow Rate

0.10-9.99 ml/hr in 0.01 ml increments
10-99.9 ml/hr in 0.1 ml increments
100-999 ml/hr in 1 ml increments

Rate Restriction by Syringe Size

50-60 ml syringe - 0.1-999 ml/hr
30-35 ml syringe - 0.1-650 ml/hr
20 ml syringe - 0.1-500 ml/hr

Occlusion Alarm Thresholds

Low, Medium, and High

Rate Accuracy

±2% of full scale plunger travel not including syringe variation