

# Alaris/Carefusion/Cardinal Medley PCA 8120

## Alaris/Carefusion/Cardinal Medley PCA 8120

### Product Specifications

#### Dimensions

4.5"W x 15.0"H x 7.5"D exclusive of security door

#### Weight

5.5 lbs.

#### Operating Principle

Positive Displacement

#### Delivery Units

mg, mcg, mL, mg/h, mcg/h, mL/h

#### Flow Rate

0.10-9.99 ml/hr in 0.01 ml increments  
10-99.9 ml/hr in 0.1 ml increments  
100-999 ml/hr in 1 ml increments

#### Rate Restriction by Syringe Size

50-60 ml syringe - 0.1-999 ml/hr  
30-35 ml syringe - 0.1-650 ml/hr  
20 ml syringe - 0.1-500 ml/hr

#### Occlusion Alarm Thresholds

Low, Medium, and High

#### Rate Accuracy

±2% of full scale plunger travel not including syringe variation