

## Alaris (Carefusion) Medley 8100

### Dimensions

3.3"W x 8.9"H x 5.5"D/p>

### Weight

2.5 lbs.

### Operating Principle

Positive Displacement

### Flow Rate Range

0.1-99.9 ml/hr in 0.1 ml increments  
100-999 ml/hr in 1 ml increments

### VTBI Range

0.1-9.99 ml in 0.01 ml increments  
10-999.9 ml in 0.1 ml increments  
1000-9999 ml in 1 ml increments

### KVO Rate

1ml/hr or selected rate, whichever is lower

### Infusion Pressure Maximum

654 mmHg

## **Rate Accuracy**

±5% at rates between 1 and 999 ml/hr

±5.5% at rates less than 1 ml/hr