

Alaris (Carefusion) Medley 8100

Dimensions

3.3"W x 8.9"H x 5.5"D/p>

Weight

2.5 lbs.

Operating Principle

Positive Displacement

Flow Rate Range

0.1-99.9 ml/hr in 0.1 ml increments
100-999 ml/hr in 1 ml increments

VTBI Range

0.1-9.99 ml in 0.01 ml increments
10-999.9 ml in 0.1 ml increments
1000-9999 ml in 1 ml increments

KVO Rate

1ml/hr or selected rate, whichever is lower

Infusion Pressure Maximum

654 mmHg

Rate Accuracy

±5% at rates between 1 and 999 ml/hr

±5.5% at rates less than 1 ml/hr