

# Alaris (Carefusion) Medley POC 8000

## **Dimensions**

6.9"W x 8.8"H x 9"D including pole clamp

## **Weight**

7.2 lbs.

## **Battery Operation**

Battery run time is a function of the number of modules attached and module activity. See operator's manual for battery operation hours.

## **Communication Data Port**

RS-232 with an RJ45 connector

## **Electronic Memory**

System configuration parameters stored in volatile memory are retained for at least 6 months by an internal backup lithium battery. Module specific parameters are stored for 8 hours when system is turned off.

## **Power Requirements**

100-240V, 50/60 Hz, 150VA Max