

Alaris (Carefusion) Medley POC 8000

Dimensions

6.9"W x 8.8"H x 9"D including pole clamp

Weight

7.2 lbs.

Battery Operation

Battery run time is a function of the number of modules attached and module activity. See operator's manual for battery operation hours.

Communication Data Port

RS-232 with an RJ45 connector

Electronic Memory

System configuration parameters stored in volatile memory are retained for at least 6 months by an internal backup lithium battery. Module specific parameters are stored for 8 hours when system is turned off.

Power Requirements

100-240V, 50/60 Hz, 150VA Max