

# Allied Autovent 3000

## Allied Autovent 3000

### Product Specifications

- Small, compact, and transport ready
- Designed to perform in harsh and toxic environments
- Assembles quickly with simplified operating instructions printed on the case
- Impact-resistant case with shock-absorbing bumpers and recessed control settings
- Rate and volume settings independently operated
- Easy-to-monitor patient valve assembly and toxic atmosphere valve are able to deliver additional breaths on demand by inspiratory effort of the patient
- Operates exclusively on standard 50 psi source gas, eliminating air entrainment and maintaining the integrity of the specified oxygen concentration
- Operates at extreme temperatures ranging from 0 F to 125 F
- Pneumatically powered, batteries not required
- Allows for PEEP attachment
- Meets ASTM, ISO and KKK-1822C standards

[Get the Quick Start Guide Here](#)