

B. Braun Infusomat® Space

Infusomat® Space, Volumetric infusion pump

A small, stackable large volume infusion pump system for all care areas that's easy to transport.

Product Description

Type of unit

Volumetric infusion pump

Dimensions

214 x 68 x 124 mm (WxHxD)

Weight

Approx. 1.4 kg

Moisture Protection

IP 22, drip protected for horizontal usage

Display

Backlit graphic display, ~40° read angle from all sides

Keypad

Backlit keys, cell phone like cursor navigation

Flow Rates

0.1-1200 ml/h

Accuracy of Set Delivery Rate

$<< \pm 0,5 \%$ mechanical accuracy

$\pm 5 \%$ according to IEC/EN 60601-2-24

Operating Temperature

+10° C ... + 40° C

+50° F +105° F

Voltage

11-16 V DC supplied by external Space Power Supply or by SpaceStation

Battery Operating Time

Minimal power consumption of the devices and new battery technology ensure long operating times, e.g. 9.4 hrs. @ 80 ml/h, 8.5 hrs. @ 100 ml/h

Lifetime

Min. 10 years under continuous duty conditions