

# **B. Braun Infusomat® Space**

## **Infusomat® Space, Volumetric infusion pump**

*A small, stackable large volume infusion pump system for all care areas that's easy to transport.*

### **Product Description**

#### **Type of unit**

Volumetric infusion pump

#### **Dimensions**

214 x 68 x 124 mm (WxHxD)

#### **Weight**

Approx. 1.4 kg

#### **Moisture Protection**

IP 22, drip protected for horizontal usage

#### **Display**

Backlit graphic display, ~40° read angle from all sides

#### **Keypad**

Backlit keys, cell phone like cursor navigation

#### **Flow Rates**

0.1-1200 ml/h

### **Accuracy of Set Delivery Rate**

$<< \pm 0,5 \%$  mechanical accuracy

$\pm 5 \%$  according to IEC/EN 60601-2-24

### **Operating Temperature**

+10° C ... + 40° C

+50° F .... +105° F

### **Voltage**

11-16 V DC supplied by external Space Power Supply or by SpaceStation

### **Battery Operating Time**

Minimal power consumption of the devices and new battery technology ensure long operating times, e.g. 9.4 hrs. @ 80 ml/h, 8.5 hrs. @ 100 ml/h

### **Lifetime**

Min. 10 years under continuous duty conditions