

# **B. Braun Outlook**

## **Mechanism**

Volumetric, Positive Pressure Displacement Reservoir

## **Dimensions**

Width - 12" Height - 4.5" Depth - 8.75"

## **Power**

120 VAC, 3 amp, 60 Hz

## **Battery**

12V Sealed lead-acid battery

## **Rate Range**

0.1-999.9 ml/hr in 0.1 ml increments  
Micro mode: 0.1-99.9 ml/hr in 0.1 ml increments

## **Volume to be Delivered Range**

0.1-9999.9 ml in 0.1 ml increments

## **Variable Occlusion Range**

## **KVO Rate**

3 ml/hr or selected rate, whichever is lower

## **Volume Delivery Accuracy**

±5%

## **Alarms**

Air-in-line

Plug in alarm

Close Roller Clamp

Container Empty

Door Open

Downstream Occlusion

Hold Time Exceeded

Low Battery Warning

System Error

## **Features**

DoseGuard™ (Outlook 100, 200)

DoseScan™ (Outlook 200)

Mode with Drug Menu  
Primary and Piggyback Modes  
Set Time and Volume  
Schedule Next Run  
Program Mode  
Variable Occlusion Pressure