

B. Braun Outlook

Mechanism

Volumetric, Positive Pressure Displacement Reservoir

Dimensions

Width - 12" Height - 4.5" Depth - 8.75"

Power

120 VAC, 3 amp, 60 Hz

Battery

12V Sealed lead-acid battery

Rate Range

0.1-999.9 ml/hr in 0.1 ml increments
Micro mode: 0.1-99.9 ml/hr in 0.1 ml increments

Volume to be Delivered Range

0.1-9999.9 ml in 0.1 ml increments

Variable Occlusion Range

KVO Rate

3 ml/hr or selected rate, whichever is lower

Volume Delivery Accuracy

±5%

Alarms

Air-in-line

Plug in alarm

Close Roller Clamp

Container Empty

Door Open

Downstream Occlusion

Hold Time Exceeded

Low Battery Warning

System Error

Features

DoseGuard™ (Outlook 100, 200)

DoseScan™ (Outlook 200)

Mode with Drug Menu
Primary and Piggyback Modes
Set Time and Volume
Schedule Next Run
Program Mode
Variable Occlusion Pressure