

B. Braun Vista Basic

Product Specifications

Mechanism

Linear peristalsis

Dimensions

5.5"W x 9.4"H x 7.8"D

Weight

6.8 lbs.

Power

100-120 VAC, 50-60 Hz, 13 Watts

Battery

7.2V, 1.2Ah NiCad rechargeable

Battery Operating Time

3.5 hours at max delivery rate

Rate Range

0.1-800 ml/hr in 0.1 ml increments with Vista™ tubing

Volume to be Infused

0.1-9999.9 ml

Occlusion Pressure

Variable Settings: Low - 8 PSI, High - 17 PSI

KVO Rate

1 ml/hr or selected rate whichever is lower

Flow Rate Accuracy

±5%

Alarms

Air-in-line

Low Battery

Plug in Now

Bag Empty

Downstream Occlusion

Door Open

Pump Stopped

KVO

Unit Malfunction

Mode/Features

Continuous and Piggyback Modes

Program Mode for variable infusions

Ramp/Taper Mode for TPN infusions

Delay A Mode for postponing start and finish