

Baxter (Bard) 300XL Mini-Infuser

Baxter (Bard) 300XL Mini-Infuser

Product Specifications

Dimensions

8.5"H x 4.5"W x 1.6"D

Weight

1 lb 12 oz with batteries

Accuracy

± -3%

Battery Life

400 hours typical, 4 hours minimum at low battery

Batteries

Four C-size alkaline cell

Flow Profile

Continuous Flat

Rate (Travel)

Xlow: 1.8 in/hr Low: 2.7 in/hr Med: 5.5 in/hr

***Consult Flow Rate Chart to determine flow rate with specific syringes**

Occlusion Force

8 lbs

***Consult Flow Rate Chart to determine occlusion pressure with specific**

syringes

Stored Volume on Occlusion

0.9 mL

Back Pressure Effect on Accuracy

None to occlusion pressure

Alarms

Normal Operation, Deliver End, Occlusion, Low Battery