

# Baxter InfusOR

## Baxter InfusOR

### Product Specifications

#### Dimensions

9.2"H x 4.5"W x 2"D

#### Weight

2lbs with batteries

#### Syringes

Monoject® or BD® 20cc or 60cc

#### Flow Rates

0-600 ml/hr

#### Flow Profile

Bolus: continuous. Infusion: pulsed intermittently

#### Max Occlusion Pressure

11 psi with 60 cc syringe; 20 psi with 20 cc syringe

#### Batteries

Four "C" Alkaline Batteries

#### Battery Life

Typical use: 150 hours

#### Accuracy

**±3%**