

Baxter InfusOR

Baxter InfusOR

Product Specifications

Dimensions

9.2"H x 4.5"W x 2"D

Weight

2lbs with batteries

Syringes

Monoject® or BD® 20cc or 60cc

Flow Rates

0-600 ml/hr

Flow Profile

Bolus: continuous. Infusion: pulsed intermittently

Max Occlusion Pressure

11 psi with 60 cc syringe; 20 psi with 20 cc syringe

Batteries

Four "C" Alkaline Batteries

Battery Life

Typical use: 150 hours

Accuracy

±3%