

## CADD Legacy Plus 6500

**Pump Resolution:**

0.050 ml per pump stroke nominal

**Dimensions:**

1.6" D x 3.8" W x 4.4" H excluding cassette

**Weight:**

13.8 oz including 2 AA batteries and empty 100 ml reservoir

**Power Source:**

2 AA alkaline batteries or AC Adapter

**Maximum Infusion Pressure:**

40.0 psi

**System Delivery Accuracy:**

±6% nominal

**Maximum Time to Occlusion Alarm:**

2.0 hours

**Delivery Rate During Priming:**

Approximately 180 ml/hr

Continuous Delivery Mode Specifications

**Reservoir Volume:**

1-9999 or Not in Use; programmable in 1 ml increments, displayed in 0.1 ml increments

**Continuous Rate:**

0.1 - 125 ml/hr; programmable in 0.1 ml/hr increments

**Given:**

0-99999.95 in 0.05 ml increments

Intermittent Delivery Mode Specifications

**Reservoir Volume:**

1-9999 or Not in Use; programmable in 1 ml increments, displayed in 0.1 ml

increments

**Dose Volume:**

0.1 - 1000 ml programmable in 0.1 ml increments

**Dose Duration:**

1 min to 24 hrs in the following increments

- 1 minute for values between 1 min and 10 min
- 5 minutes for values above 10 min

**Dose Cycle:**

10 min to 96 hrs in 5 min increments

**KVO Rate:**

0 to 125 ml/hr in 0.1 ml/hr increments

**Dose Starts In:**

Immediate or 1 min to 96 hrs in the following increments

- 1 min for values between 0 min and 10 min
- 5 min for values between 10 min and 96 hrs

**Given:**

0-99999.95 in 0.05 ml increments

**Dose Remaining:**

1 minute increments

**Pump Alarms:**

Low battery power

Depleted battery power

Battery dislodged

Pump stopped

Pump fault

Low reservoir volume

High delivery pressure

Air in line

Disposable not attached when run attempted

Motor locked

Upstream occlusion

Reservoir volume empty

Program incomplete

Key stuck

Disposable detached