

CareFusion LTV 950

CareFusion LTV 950

Product Specifications

Dimensions

Height 3 inches (8 cm) 3.25 inches (8.4 cm) Width 10 inches (25 cm) 10.5 inches (27 cm) Depth 12 inches (30 cm) 13.5 inches (38 cm)

Weight

Bare 12.85 lbs (5.83 kg) Loaded 14.45 lbs (6.5 kg)

Power indicators

External power: Green(Full Power), Amber(Low Power) Charge status: Flashing Amber(Pretest), Green(Charge Complete), Amber(Charging), Red(Charger or Battery Fault) Battery level: Green(Acceptable), Amber(Low), Red(Critically Low)

Monitors and Indications

Peak inspiratory pressure: 0 to 120 cmH₂O Mean airway pressure: 0 to 99 cmH₂O PEEP: 0 to 99 cmH₂O Total breath rate: 0 to 250 bpm Airway pressure display: -10 to 108 cmH₂O Inhale/Exhale ratio: 99:1 to 1:99 Calculated peak flow: 10 to 100 lpm Patient effort: Green LED

Alarms

Disconnect/sense line Low and lost external power Low and empty internal battery Apnea interval 10 to 60 seconds High pressure limit 5 to 100 cmH₂O Low peak pressure Off, 1 to 60 cmH₂O Alarm volume 60 to 80 dBA at 1 meter Alarm silence/reset 60 seconds Ventilator inoperative Red LED, audible alarm

Power Cord

Biomed replaceable 2.75 m (9 ft)

Extended Features

NPPV (Non-Invasive Positive Pressure Ventilation) mode Apnea back-up ventilation

The LTV-Series of Ventilators all look, act, and, feel alike with the exception of which controls and features are installed in each model. The following list is unique to this model.

- Pressure control
- Volume control
- Pressure support
- Flow-Triggering
- Leak Compensation
- Multiple power source
- Non-invasive ventilation
- **Low pressure O2 inlet**

[Get the User Manual Here](#)