

Datex/Ohmeda: Aestive/5

Superior ventilation:

Volume mode, pressure mode, electronic PEEP
Tidal volume compensation
One motion from mechanical to manual mode
Two key presses to total standby: end case
Cardiac bypass case mode

Open systems architecture:

Lower overall height
User configurable drawers/shelving

Innovative patient breathing system:

Eight machine hoses/cables integrated
“No tools” disassembly of components
Autoclavable and latex-free
Responsive location of common gas outlet

Improved low flow/reduced life cycle costs:

Fresh gas flow compensation-automatically
Smooth, faster acting fresh gas flow control
Minimum O2 flow of 50 mL
Dual air flow tube for low flow