

Datex/Ohmeda: Aestive/5

Superior ventilation:

Volume mode, pressure mode, electronic PEEP

Tidal volume compensation

One motion from mechanical to manual mode

Two key presses to total standby: end case

Cardiac bypass case mode

Open systems architecture:

Lower overall height

User configurable drawers/shelving

Innovative patient breathing system:

Eight machine hoses/cables integrated

“No tools” disassembly of components

Autoclavable and latex-free

Responsive location of common gas outlet

Improved low flow/reduced life cycle costs:

Fresh gas flow compensation-automatically

Smooth, faster acting fresh gas flow control

Minimum O₂ flow of 50 mL

Dual air flow tube for low flow