

## Dräger Babylog 8000 plus

### Dräger Babylog 8000 plus

#### Features:

- Flow and Volume Monitoring at the Y-piece
- Fully integrated airway monitoring and graphics system which displays waveforms for both pressure and flow
- Flow-triggered ventilation, measured at the Y-piece offers more predictability and control
- Digital valve array of the Babylog allows adjustable flow patterns.
- Changes in ventilator settings can be seen immediately in real time display of pressure and flow wave forms and from measured volume values
- Air trapping can be seen at a glance.