

# Dräger Evita Infinity V500

## Dräger Evita Infinity V500

### Benefits

#### Tools for your ventilation therapy I

- Lung diagnostic tools like the Low Flow maneuver (inflection points)
- Breath by breath recruitment trends (e.g. PEEP, EIP, VT, Cdyn)
- Recruitment tools (e.g. Inspiration Hold, QuickSet, PressureLink)
- PC-APRV with AutoRelease
- Volumetric CO<sub>2</sub>-Monitoring (VCO<sub>2</sub>, VTCO<sub>2</sub>, Slope Phase 3, V<sub>d</sub>s/V<sub>T</sub>e)
- Weaning parameter (e.g. RSBi, P0.1, NIF)
- Automatic weaning with SmartCare/PS®

#### Tools for your ventilation therapy II

- Variable Pressure Support, Proportional Pressure Support, Automatic Tube Compensation
- Graphical representation of airway resistance and lung compliance with Smart Pulmonary View
- “Room-to-breathe” concept (AutoFlow®, BIPAP, VG)
- Applicable for neonatal ventilation (smallest tidal volume: 2 ml)

#### Functions to support your workflow

- Non-invasive ventilation in all modes and all patient categories
- O<sub>2</sub>-Therapy allows flow constant oxygen application
- Flexible screen configuration: 6 different views for each patient related to the individual therapy

- Full record of all patient data, alarms and trends
- Data export via USB interface
- Context-sensitive help function and online Instruction For Use for ventilation modes, alarms and device functions
- RFID functionality e.g. to monitor exchange intervals and to transport ventilation settings