

# Dräger Evita XL

## Dräger Evita XL

The Dräger Evita® XL delivers a comprehensive array of both controlled and assisted ventilation modes. It combines a broad range of performance capabilities with many advanced application functionalities, providing outstanding ventilation therapy and easy-to-use operation encompassing the entire scope of the ventilation process. The Evita® XL helps caregivers find optimal ventilation settings, supports the weaning process, and improves caregiver workflow.

## Benefits

### Non-Invasive NIV Plus

Features an advanced, dynamic leak compensation system that provides adaptive responsiveness and reliable tidal volume delivery.

### Open breathing system with standard AutoFlow®

Standard AutoFlow® and BIPAP1) /PCV+, makes spontaneous breathing possible in any ventilation mode. Eases change to a stable breathing patient.

### Help prevent complications

The Lung Protection Package (LPP) provides you with a low flow P/V manoeuvre to identify optimal-PEEP and pressure settings for protective ventilation while special recruitment trends monitor the process. Integrated CO2 monitoring helps verify correct intubation and metabolic stability.

### AutoFlow™

Automatic adaptation of inspiratory flow in volume controlled modes .

## Specifications

### Dimensions

530 x 315 x 450 mm / 20.9 x 12.4 x 17.7 inches (without trolley)

### Weight

Approx. 50 kg / 64 lbs

[View the technical data here.](#)