

Drager: Narkomed 2B

Features:

ORM standard; ORMC optional.

Ventilator options: Pneumatic or electronic, ascending or descending bellows.

Oxygen monitor and CO2 absorber standard.

AVE Electronic ventilator.

Absorber and scavenger

Oxygen analyzer

Breathing pressure monitor

Respiratory volume

Centralized alarm and data

Oxygen ratio controller

Pressure limit controller

Drager 19.1 vaporizer available

Specifications:

Maximum Height: 68" (173cm)

Maximum Width: 40" (102cm)

Maximum Depth: 25" (63.5cm)

Table-top Height: 34" (86cm)

Table-top Area: 271 in. sq. (1,749cm²)

Shelf Area: 200 in. sq. (1,290cm²)

Drawer Dimensions:

Large: 14"x21"x7" (36x52x18cm)

Small: 14"x17"x3" (36x43x8cm)

Nominal Weight: 325 lb. (147kg)