

# Drager: Narkomed 4

## Features:

19.1 vaporizers  
AVE ventilator  
NIB  
SPO2  
CO2  
Built in monitor

## Options:

Air system  
PEEP  
Boom Arm

## Specifications:

Weight: 500lb (226.8kg)  
Height: 64 3/4 in (165.5 cm)  
Width: 31 3/4 in (80.6 cm)  
Depth: 25 1/8 in (63.8 cm)  
Shelf area: 200 in sq (1,290 cm sq)  
Power Supply: 117 volts  
Optional Power Supply: 220/240 volts  
Backup Battery: Charging time: <12 hours Reserve power time: >25 minutes