

Drager: Narkomed 4

Features:

19.1 vaporizers
AVE ventilator
NIB
SPO2
CO2
Built in monitor

Options:

Air system
PEEP
Boom Arm

Specifications:

Weight: 500lb (226.8kg)
Height: 64 3/4 in (165.5 cm)
Width: 31 3/4 in (80.6 cm)
Depth: 25 1/8 in (63.8 cm)
Shelf area: 200 in sq (1,290 cm sq)
Power Supply: 117 volts
Optional Power Supply: 220/240 volts
Backup Battery: Charging time: <12 hours Reserve power time: >25 minutes