

# GE: Aestiva 5/7100 Ventilator

## Features:

Low overall height

Exhaled volume, airway pressure and inspired oxygen monitoring capabilities

Direct access to ventilator parameters

Pressure bar graph for visual reference on a breath-by-breath basis

Mechanical bag/vent switch turns ventilator on/off

Smart alarms direct user to specific problems

Fresh gas flow compensation- automatically

Autoclavable and latex-free

Integrated machine hoses/capabilities

## Dimensions:

Height: 134.5 cm/52.9 in

Width: 72 cm/28.3 in

Depth: 73 cm/28.7 in

Weight: Approximately 108 kg/238 lb

## Top shelf:

Weight limit: 34 kg/75 lb

Width: 66 cm/26 in

Depth: 40 cm/15.75 in