

Gomco/Allied Medical: Constant/Intermittent Pump Model 6036

Gomco/Allied Medical: Constant/Intermittent Pump Model 6036

Product Specifications

Vacuum Range

250 mm Hg (constant) / 200 mm Hg (intermittent)

Flow Rate

2.5 lpm (constant) / 3.4 (intermittent)

Collection Bottle

Provided with two glass collection bottle systems, one for each mode selection (one 2800 ml and 600 ml)

Shipping Weight

78 lb

Dimensions

33"H x 23"W x 17"D

Compressor

Diaphragm

Motor

1/8 hp Continuous Duty

Features

- Mode Vacuum Range Adjustable On/Off Cycle Time
- Constant 0-250 mm HG On — 5-90 seconds
- **Intermittent 0-200 mm Hg Off- 23 seconds**