

# Kendall 7325

## Kendall 7325

### Product Specifications

#### Dimensions

6"H x 10.75"W x 4.5"D

#### Weight

9.3 pounds

#### Power Requirements

115 VAC, 50/60 Hz

#### Compression Type

Sequential, Gradient Pressure

#### Compression Cycle

11 seconds compression, decompression time based upon Vascular Refill Detection Measurement

#### Set Pressure

45 mmHg

#### Sleeve Cooling

Yes

#### Cycle Monitor

Yes

#### Bed Hook

**Yes**

**Alarms**

**Low Pressure, High Pressure, Internal Electronics Malfunction**