

# Kinetec Spectra Knee CPM

An unbeatable combination of good looks and mechanical innovation sets Spectra apart as the Knee CPM of choice. The entire spectrum of needs for knee patients can now be addressed with a carriage that accommodates a full range of patients, without any additional accessories. Lightweight and exceptionally easy to carry, Spectra offers a quick set-up and practical operational modes that enhance patient comfort, compliance and results.

## Features

- Long-life DC motor and screw-drive mechanism
- Anatomically correct alignment
- Wide carriage (10") accommodates all patients
- Hyperextension (-10°) to full Flexion (120°)
- Internally contained medical grade transformer
- Exclusive Load Reverse ensures patient protection
- Lightweight: only 24 pounds
- Sleek design incorporates dual carry handles
- Comfortable, easy-to-apply patient pad kit included
- Patient lock-out prevents accidental program changes
- Waterproof digital hand control that snaps conveniently into the foot support
- Pause setting for up to 15 minutes of rest

## Product Specifications

- Length: 37" (94cm)
- Full Leg Size: 28" to 41" (71 to 104 cm)
- Tibia: 15" to 23" (38 to 58 cm)
- Femur: 13" to 18" (33 to 46 cm)