

# Medfusion 3500

## Medfusion 3500

### Product Specifications

#### Dimensions

10.5"W x 5.1"H x 5.5"D

#### Weight

3.8 lbs without pole clamp

#### Accuracy

± 2% excluding syringe variations

#### Flow Rate Range

0.01 ml/hr to 1130 ml/hr depending on syringe size

#### Delivery Modes

ml/hr, body weight, mass, volume over time, customer dilution, and intermittent

#### Delivery Options

Loading Dose, Bolus Dose, Standby, Volume Limit, and KVO Rate

#### Power

100 to 240 VAC, 50/60 Hz, 30VA

Internal rechargeable battery

Optional external DC Input, 6-14 VDC, 7 Watts

#### Safety Class II - Type CF

#### Internal Battery Capacity

**10 hours at 5.0 ml/hr with 60 ml syringe**

**User Defined Delivery Libraries**

**Four Libraries with 16 profiles per library**

**MR Tested**

**May be used in MRI Rooms with static magnetic fields up to 150 gauss**

**Options**

**PharmGuard™ Error Reduction Software for programmable limits**