

Medfusion 3500

Medfusion 3500

Product Specifications

Dimensions

10.5"W x 5.1"H x 5.5"D

Weight

3.8 lbs without pole clamp

Accuracy

± 2% excluding syringe variations

Flow Rate Range

0.01 ml/hr to 1130 ml/hr depending on syringe size

Delivery Modes

ml/hr, body weight, mass, volume over time, customer dilution, and intermittent

Delivery Options

Loading Dose, Bolus Dose, Standby, Volume Limit, and KVO Rate

Power

100 to 240 VAC, 50/60 Hz, 30VA

Internal rechargeable battery

Optional external DC Input, 6-14 VDC, 7 Watts

Safety Class II - Type CF

Internal Battery Capacity



10 hours at 5.0 ml/hr with 60 ml syringe

User Defined Delivery Libraries

Four Libraries with 16 profiles per library

MR Tested

May be used in MRI Rooms with static magnetic fields up to 150 gauss

Options

PharmGuard™ Error Reduction Software for programmable limits