

Newport HT50

Newport HT50

The Newport HT-50 All-Purpose ventilator offers comprehensive, versatile features to meet the demands of almost any healthcare setting. The HT-50 is suitable for both pediatric and adult patients in hospital, extended care, home care or transport use.

Features

- Self-contained Volume or Pressure Ventilator
- Sophisticated Modalities such as A/CMV, SIMV, SPONT
- Volume or Pressure Control/Pressure Support
- PEEP/CPAP
- Comprehensive Alarm functions
- Built-in Battery (up to 10 hr. operation) with integrated alarm system
- Self-Switching Versatile Power Supply
- Patented self-contained air compressor system. External compressed air is not needed.
- Brightly lit LED displays identify essential functions

[Get the User Manual Here](#)