

# Physio Control Lifepak 12 Defibrillator

The Lifepak12 Biphasic defibrillator/monitor series provides therapeutic and diagnostic functions in a single small device designed for both out-of-hospital and hospital users. Innovative platform design provides ease of use and the flexibility to add new features and enhancements at a future date.

Fast, safe and effective defibrillation and pacing therapy is delivered through single function keys. The Selector knob allows easy input of events and patient demographics and facilitates navigation through menu options.

In AED mode, the LIFEPAK 12 defibrillator/monitor series utilizes the Physio-Control field proven Shock Advisory System to advise operator when it detects a shockable rhythm. In manual mode, the LIFEPAK 12 defibrillator/monitor series features simple 1-2-3 operation of synchronized cardioversion.

Large LCD display allows 1, 2 or 3 ECG channels to be viewed simultaneously, with up to 8 seconds of cascading ECG.

Pulse oximetry is available with Nellcor C-LOCK ECG synchronisation for timing pulse oximetry measurements with ECG signal.

The optional LifePak 12 AC Power Adapter (sold separately) provides line power as well as battery charging capability.