

Precision Aspirator

Medical:

EasyVac

Precision Medical: EasyVac Aspirator

Product Specifications

Dimensions

Length 13 1/4" (33.7 cm), Width 8" (20.3 cm), Height 9 3/4" (24.8 cm)

Weight

12.6 lbs. (5.7 kg)

Shipping Weight

14.5 lbs. (6.6 kg)

Motor

120 VAC, 60Hz, 5 Amps

Vacuum Range

0-26 inHg

Horsepower

1/5

Operating Temperature Range

50°F to 104°F (10°C to 40°C)

Storage Conditions

Temperature: -4°F to 104°F (-20°C to 40°C) Humidity: Max 95% non-condensing