

# Respironics Trilogy 100

## Respironics Trilogy 100

### Simplicity

Light, versatile, and easy to use, the Trilogy 100 offers clinicians and patients one of life's greatest qualities – simplicity. Whether at home or alternative-care sites, the Trilogy 100 portable volume and pressure support ventilator with BiPAP technology is made for adult and pediatric (>5kg) use.

### Two Circuit Options

Interchangeable passive and active exhalation porting blocks give clinicians the flexibility to choose the best circuit and interface available for patients. Its internal Peep, combined with a single limb circuit, creates an uncomplicated application.

#### Passive Exhalation Porting Block

- Uses a single BiPAP circuit with Whisper Swivel II or a mask with an integrated leak
- **Is compatible with pressure and volume ventilation modes**

#### Passive Exhalation Porting Block

- Uses a single limb circuit with a proximal pressure line and active exhalation valve.
- **Is compatible with pressure and volume ventilation modes**

### Portability

- Trilogy 100 was designed with maximum portability in mind
- Only weighs 5kg or 11lbs

- Has 3 hour internal battery and 3 hour detachable battery
- **Able to connect to an external lead acid battery**

## Data Management

As with many of our products, The Trilogy 100 comes equipped with an SD Data card to allow clinical evaluation of use and patient progress.

Get the Quick Start Guide Here