

# Respironics Trilogy 100

## Respironics Trilogy 100

### Simplicity

Light, versatile, and easy to use, the Trilogy 100 offers clinicians and patients one of life's greatest qualities – simplicity. Whether at home or alternative-care sites, the Trilogy 100 portable volume and pressure support ventilator with BiPAP technology is made for adult and pediatric (>5kg) use.

### Two Circuit Options

Interchangeable passive and active exhalation porting blocks give clinicians the flexibility to choose the best circuit and interface available for patients. Its internal Peep, combined with a single limb circuit, creates an uncomplicated application.

#### Passive Exhalation Porting Block

- Uses a single BiPAP circuit with Whisper Swivel II or a mask with an integrated leak
- **Is compatible with pressure and volume ventilation modes**

#### Passive Exhalation Porting Block

- Uses a single limb circuit with a proximal pressure line and active exhalation valve.
- **Is compatible with pressure and volume ventilation modes**

### Portability

- Trilogy 100 was designed with maximum portability in mind
- Only weighs 5kg or 11lbs

- Has 3 hour internal battery and 3 hour detachable battery
- **Able to connect to an external lead acid battery**

## Data Management

As with many of our products, The Trilogy 100 comes equipped with an SD Data card to allow clinical evaluation of use and patient progress.

[Get the Quick Start Guide Here](#)