

# Respironics Universal Rolling Stand

Respironics Universal Rolling Stand for the V60 Ventilator or BiPAP Vision (Pictured with BiPAP Vision)

This is a Respironics Universal Rolling Stand that can be used for the \*V60 Ventilator or BiPAP Vision. There is a different mounting plate for each. Please specify which one you will need when ordering.

*\*Compatible with V60 ventilators with serial numbers beginning with "1": all serial numbers below 100111553; V60 ventilators with serial numbers beginning with "2": all serial numbers below 201007805.*

## **Product Specifications**

- Built-in humidifier bracket
- O2 hose mount
- Power cord mount
- Accessory basket
- Large castors with brakes

Drop ships directly from manufacturer within 1-7 business days.