

Schuco: Vac 330 Aspirator Suction Unit

Schuco: Vac 330 Aspirator Suction Unit

The Schuco-Vac 330 is a durable medical suction unit which provides variable controlled vacuum regulation between 0-22 Hg and a suction flow rate of 28 lpm. Features vibration-free, easy to read regulator gauge and a hydrophobic/bacteria filter to prevent pump flooding and contamination. Accepts 800cc & 1200cc canisters. 1/10HP motor, 115V 60Hz.

Product Specifications

Vacuum Regulation

0-22 Hg

Suction Flow Rate

28 lpm

Collection Canister

Accepts 800cc and 1200cc

Motor

1/10HP, 115V 60Hz