

# Schuco: Vac 330 Aspirator Suction Unit

## **Schuco: Vac 330 Aspirator Suction Unit**

The Schuco-Vac 330 is a durable medical suction unit which provides variable controlled vacuum regulation between 0-22 Hg and a suction flow rate of 28 lpm. Features vibration-free, easy to read regulator gauge and a hydrophobic/bacteria filter to prevent pump flooding and contamination. Accepts 800cc & 1200cc canisters. 1/10HP motor, 115V 60Hz.

## **Product Specifications**

### **Vacuum Regulation**

0-22 Hg

### **Suction Flow Rate**

28 lpm

### **Collection Canister**

Accepts 800cc and 1200cc

### **Motor**

1/10HP, 115V 60Hz