

# Smiths Medical CADD 100 mL Medication Cassette Reservoir 21-7302-24

## Smiths Medical CADD 100 mL Medication Cassette Reservoir

CADD Medication Cassette Reservoir, 100 mL (Clear), with Flow Stop, Clamp,  
Female Leur, Latex Free

All medication cassette reservoirs with flow stop free-flow protection meet the JCAHO National Patient Safety Goals for prevention of free-flow of medication. Flow stop free-flow protection is an integral set-based free-flow protection device that has been added to medication cassette reservoirs.

Order quantity: 12/box

Shelf life: 2 years

Priming volume: 0.2 mL

Rates up to: 125 mL/hr

Tubing length: 8"

Plasticizer: DEHP

Please see the *Instructions for Use* for a complete listing of the indications, contraindications, warnings and precautions.

### PRODUCT BENEFITS

- Small
- Lightweight
- Portable