

# Smiths Medical CADD 250 mL Medication Reservoir

## Smiths Medical CADD 250 mL Medication Reservoir

CADD Medication Cassette Reservoir, 250 mL (Clear), with Flow Stop, Clamp, Female Leur, 8" tubing

The CADD™ 250 mL medication cassette reservoirs are designed to provide safer, more secure, ambulatory medication delivery for infusion therapies requiring volumes up to 250 mL. CADD™ medication cassette reservoirs promote patient ambulation and minimize the impact of infusion therapy on a patient's normal lifestyle by decreasing the length of post-operative hospitalizations and time away from home.

Order quantity: 12/box

Shelf life: 2 years

Priming volume: 0.2 mL

Rates up to: 125 mL/hr

Tubing length: 8"

Plasticizer: TOTM

Please see the *Instructions for Use* for a complete listing of the indications, contraindications, warnings and precautions.

### PRODUCT BENEFITS

- Small
- Lightweight
- Portable
- Lockbox alternative