

Smiths Medical CADD 250 mL Medication Reservoir

Smiths Medical CADD 250 mL Medication Reservoir

CADD Medication Cassette Reservoir, 250 mL (Clear), with Flow Stop, Clamp, Female Leur, 8" tubing

The CADD™ 250 mL medication cassette reservoirs are designed to provide safer, more secure, ambulatory medication delivery for infusion therapies requiring volumes up to 250 mL. CADD™ medication cassette reservoirs promote patient ambulation and minimize the impact of infusion therapy on a patient's normal lifestyle by decreasing the length of post-operative hospitalizations and time away from home.

Order quantity: 12/box

Shelf life: 2 years

Priming volume: 0.2 mL

Rates up to: 125 mL/hr

Tubing length: 8"

Plasticizer: TOTM

Please see the *Instructions for Use* for a complete listing of the indications, contraindications, warnings and precautions.

PRODUCT BENEFITS

- Small
- Lightweight
- Portable
- Lockbox alternative